

# MND TAS equipment library loan scheme – Equipment list (updated May 2021)

## Mainstream Eye Gaze System

- Suitable for use as communication device and only simple internet use for people requiring eye gaze accessible technology.
- Tobii Eye Tracker - Tracker 4C
- Set up with Surface Pro tablet and includes desktop mounting arm.
- software/apps for communication and internet use

## High End Eye Gaze System

- Suitable for use as communication device and internet use for people requiring eye gaze accessible technology.
- Tobii Dynavox PC Eye Mini featuring classic gaze interaction software with mouse emulation and gaze selection
- Set up with Surface Pro tablet and includes desktop mounting arm
- Tobii Dynavox Gaze Interaction Software – Other software/apps for communication and internet use

## Switch Accessible Communication System

- Suitable for use as communication device and internet use for people requiring switch accessible technology.
- JoyCable lets you attach one or two 3.5mm switches to the tablet via a USB 2.0 connection.
- Set up with Surface Pro tablet and includes desktop mounting arm.

## Head Tracking Communication System

- Suitable for use as communication device and internet use for people requiring head tracking accessible technology
- Tracker Pro mouse replacement device which allows people to access the tablet via head movements, connected via USB port.
- Set up with Surface Pro tablet and includes desktop mounting arm

## Wheelchair, floor, desk mounting systems and mounting plates

- Suitable for use to mount switches, Surface Pro tablets, and iPads?

## iPads

- Suitable for use with clients for communication and internet access who can use touch screen controls. Not currently switch accessible.
- 10.2" 32GB Wi-Fi tablets
- Predictable communication app loaded.
- iAdapter 5A for iPad Air x 1 – Available as add on for use to enhance use of the iPad as a communication device by making it more durable for portability, amplifying the sound, and improving audio quality.

## Lightwriters

- Lightwriter SL35M
- Portable communication device for use with clients who can use Lightwriter keyboard. Features include: dual display with text displayed facing the user and the communication partner; DECTalk speech synthesis with English male, female and child voices available; memory of 36 commonly used phrases available; word prediction; and tremor control.
- Scanning module Lightwriter SM
- Scanning module available to convert standard Lightwriter to one with scanning access for clients who cannot use the keyboard. Can be used with 1 or 2 switches.

## Environmental Control Equipment

- Google Home Mini – A voice-controlled speaker and environmental control unit that can be used to play music, control smart home gadgets, add items to a shopping list, create calendar appointments, play video on a Chromecast enabled screen, etc.
- Amazon Echo Show 8 – Voice-controlled speaker and environmental control unit with 8" screen. Enables video calls to be made when text to speech software is in use.
- Broadlink RM Pro – Environmental control unit used through tablet or phone application. Allows control of infrared/RF appliances through smartphone.

## Other Equipment

- Voice amplifiers
- Light articulated switch mount, multi-clamp, and jellybean adaptor- Available to mount switches to multiple different surfaces using the multi-clamp.
- Portable hard cases for Surface Pro tablets
- Laptop soft bags

## Personal Alarms

- 4G Live Life personal falls detector alarms – to be set up with clients' emergency contacts.

To progress the order for the LiveLife personal alarm with falls detection, contact details are required to program the alarm with contacts.

The personal contact can use multiple number entries (e.g. Jane Smith work and Jane Smith mobile). It is a maximum of six numbers per device. There is no obligation to use all six contacts.

Discuss the preferred contacts with family/neighbours and send the completed table back to ILCT for processing.

Name:			
Home address:			
Home phone number:		Colour	Black/White
Contact 1 name:		Contact 1 phone number:	
Contact 2 name:		Contact 2 phone number:	
Contact 3 name:		Contact 3 phone number:	
Contact 4 name:		Contact 4 phone number:	
Contact 5 name:		Contact 5 phone number:	
Contact 6 name:		Contact 6 phone number:	